



1. Intraarticular shaver system with cutter blade used in resection of anterior flap in complex tear of lateral meniscus (under television observation)



2. Knife blade inserted through cannula, cutting (from anterior to posterior aspect to trim anterior horn)



3. Basket forceps removes segment of torn posterior horn after section by scissors



4. Final cleanup of meniscus with trimmer



5. Immediate postoperative appearance after release of tourniquet



6. Three-month postoperative arthroscopic examination. Patient asymptomatic. Note fibrous tissue response at margin of meniscus